

instructions for use
EKG-Patch for long-term ECG monitoring



Rev.1

General Information / Safety Instructions

The instructions for use must be read carefully before the first start-up and kept for future reference. Non-compliance or improper use releases the manufacturer from liability. Before each use, the medical device and its accessories must be checked for function and intactness. FAROMED GmbH Medizintechnik expressly warns against changing the medical device and its accessories. Any change leads to the exclusion of liability. The CE mark expires with the use of non-approved spare parts. The product may only be used by medically trained personnel!

Warnings

Reusing Zeus without permission can lead to cross-infection as cleaning residue or other contaminants may adhere to it. Make sure that the product is intact, damaged devices cannot be used. Protect the device against electrical currents, fields and radiation as well as high temperatures. The battery poses a risk of explosion. Only combinations approved by the manufacturer are allowed. The device can be damaged if an external defibrillator is used. This EKG patch is not waterproof and must be protected from moisture penetration.

Handling

Protect the product from mechanical damage! Do not throw! Do not use any force! Apply the patch the right way round. For more information on handling, see the body sensor instructions for use.

Disposal

After expiration of usage the product needs to dispose in accordance with your local laws, especially the battery. For the European Union is the following directive valid: According to the European directive 2012/19/EU on electrical and electronic equipment (WEEE) is this product not allowed to get disposed in your household waste. Prepare the product and its accessories for recycling or separate collection according to the directive. The directive is only valid for not contaminated products.

Symbols and inscriptions

	CE labelling The device complies with the requirements of the European Council Directive 93/42 / EEC for medical devices.
	Information about the manufacturer Shows the manufacturer of the medical device.
	hazardous waste Must not be disposed of with normal household waste; Must be disposed of in an environmentally friendly manner in accordance with local regulations.
	do not reuse Refers to a medical device that is intended for single use or for use on a single patient during a single treatment (e.g. EKG patch).
	Attention - note the accompanying documents! Refer to the accompanying documents before use and if anything is unclear.

	unsterile To indicate that the device that is normally provided sterile in the same or similar packaging has not been sterilized.
	Humidity, limitation Describes the moisture range to which the medical device can be safely exposed.
	temperature limitation The upper and lower temperature limit values must be indicated next to the upper and lower horizontal lines.
	Keep dry
	Do not use if the packaging is damaged Indicates a medical device that should not be used if the packaging is damaged or opened.

Intended use

The EKG patch ZEUS is only intended for attaching the EKG recorder ARES from AthenaDiaX to the patient.

Note:

Replace the monitoring extension (ECG electrodes) every 24 hours.

Location, duration, frequency of use

This patch is not reusable and is worn for one, two or up to seven days, depending on the variant.

Contraindication

See instructions for use body sensor.

Technical data

biocompatibility of the adhesive film	DIN EN10993-1; DIN EN10993-10
distribution	AthenaDiaX GmbH, Heinrich-Mann-Allee 3b 14473 Potsdam
size	57,4 x 57,4 x 10 mm
weight	Ca. 24 g
batteries	Lithium-Manganese-Dioxid; 3,0 V; 450 mAh - 700mAh
durability	2 years, at a temperature of 21°C and relative humidity of 50% in the original packaging
classification	Class I acc. Annex IX, RL 93/42 EWG
Item number	Zeus 1D: REF 01-102-1D Zeus 2D: REF 01-102-2D Zeus 7D: REF 01-102-7D

	Observe the operating instructions Indicates the need for the user to refer to the instructions for use.
	batch code Displays the manufacturer's batch designation so that the batch or lot can be identified.
	item number Displays the manufacturer's order number so that the medical device can be identified.
	serial number Displays the manufacturer's serial number so that a specific medical device can be identified.
	application part: Type BF
	contains or presence of natural rubber latex Indicates the presence of natural rubber or dry natural rubber latex as a construction material in the medical device or the packaging of a medical device.

Support

For questions and problems with the body sensor please contact:

AthenaDiaX GmbH
E-Mail: Info@athenadiax.de
www.athenadiax.de
+49 30/ 692 067 950



Faromed GmbH, Saalmanstrasse 9,
13403 Berlin, Deutschland

